



#### **MODELS** 2-WALL GRAB RAIL & BACKREST

RBA4150-229 90° Rail, Polyurethane Backrest, Stainless Steel Grab Rail I RH RBA4150-230 90° Rail, Polyurethane Backrest, Stainless Steel Grab Rail I LH RBA4150-239 40° Rail, Polyurethane Backrest, Stainless Steel Grab Rail I RH RBA4150-240 40° Rail, Polyurethane Backrest, Stainless Steel Grab Rail I LH



\*RBA4150-229 Shown



As improvements in the design and performance of RBA products are continuous, specifications may be subject to change without notice. The illustrations and descriptions herein are applicable to production as of the date of this Installation Instructions Sheet. Revised 12/22 © 2022 by RBA Group II/MODEL RBA4150-2 series/1222



#### **Installation Instructions**


- 1. Mount Backrest concealed wall fixing plate in a location as per Fig. A
- 2. Assemble Polyurethane Backrest to rail with hardware provided.
- 3. Based on location of toilet seat hinge point, determine desired length of the Backrest to achieve 95°-100° angle on seat when raised as shown in Fig B.
- 4. Cut stainless steel bar of Backrest to suit the required angle of the toilet seat when open.
- 5. Loosely fit mounting brackets for adjoining rails over the ends of the Backrest rail. The FxF connector [Fig. C] should be on the rail near the side wall and the MxF connector [Fig. D] on the opposite side.
- 6. Slide Backrest mounting brackets with dress rings onto the stainless steel bar of the Backrest. Tighten the mounting bracket cover with the tension bolts to a torque of 80Nm.
- 7. Set the bracket cover over concealed mounting bracket and secure with locking bolts on the bottom of the mounting bracket cover.
- 8. Set toilet seat upright to confirm correct position.
- 9. Slide the side wall rail into the bracket located on the Grab Rail.
- 10. Locate the side wall rail in the correct location in accordance with the drawing & AS1428.1
- 11. Remove cover flange and mount unit to wall. Fixings by others.
- 12. Tighten allen head grub screws on FxF connector.
- 13. Slide the rear wall rail into the MxF connector.
- 14. Remove cover flange and mount unit to wall. Fixings by others.
- 15. Tighten allen head grub screws on the MxF connector.
- 16. Replace covers over the rail mounting flanges ensure raised edges of the flange are lined up with the dimples of the bracket cover and press into place.


### **Notes**

- 1. Measure site to find correct location for Grab Rail & Backrest ensuring it is installed to meet AS1428.1 and local building
- 2. Ensure wall construction is able to support the unit and loads it is subject to. Grab Rails & Backrests are no stronger than the anchors & walls to which they are attached.
- 3. Wall anchors to be supplied by others.
- 4. Rough-in should be done with each unit. Construction is subject to manufacturers tolerances of 10mm.



# Rough in





#### **Rail Connector**




Fig. C: Male x Female (MxF) Connector




Fig. D: Female x Female (FxF) Connector

# **Cleaning & Maintenance Guide**

- 1. To remove water spots or rust spots, stainless steel cleaner/polish on a cloth is recommended.
- 2. If there are stubborn spots or if you wish to treat a scratch, synthetic abrasive general purpose pads such as scotch brite are recommended.
- 3. Apply stainless steel cleaner/polish to the synthetic abrasive pads and carefully rub the panel with the grain.
- 4. DO NOT use harsh chemicals, abrasive or petroleum based cleaners. Use of these will void the RBA warranty. DO NOT use abrasives on powder coated units.
- 5. Stainless steel should be kept clean at all times. If a coating of stainless steel cleaner/polish is maintained, stainless steel surfaces will retain their new, clean, polished appearance indefinitely. Use clean mild soapy water for powder coated units.






Fig. E: RBA4150-229 90° Right Hand 2 Wall grab Rail With Backrest








Fig. F RBA4150-230 90° Left Hand 2 Wall grab Rail With Backrest









Fig. G RBA4150-239 40° Right Hand 2 Wall grab Rail With Backrest










Fig. H RBA4150-240 40° LeftHand 2 Wall Grab Rail with Backrest

**Note:** This product should be installed, by suitably qualified persons, in a fit for purpose application, to suitable materials, using suitable fixings and comply with any relevant codes. It should be inspected periodically for signs of wear and tear that may affect performance or safety.

Dimensions are subject to manufacturer's tolerance of +/-10mm. Rough-in should be completed with each fixture.

**Important:** Installation Instructions are subject to change without notice. Please visit **our websites** for latest revision.

