

## Model RBA4007-300-Series



\*RBA4007-312 Shown

### Materials & Features

Suitable for Bariatric assistance. Rail can be folded up against the wall in a simple one-hand operation, to allow easy access to the WC. Auto-Hold device retains the rail in place when returned to the upright position. Includes lockable support leg with height adjustable slip resistant footing. Nurse call and/or toilet roll holder are welded to the main body of the rail and cannot be moved.

### Technical Detail

Bar:	Type-304 Stainless Steel, 38mm Dia.
Finish:	Gloss Grey Powder Coat
Mounting plate:	5mm Stainless Steel
SWL:	500kgs downwards 160kgs lateral
Fixings:	By others, RBA recommends M12 or equivalent screw to suit install location.

### Installation

Warning: Grab Rails are no stronger than the anchors or walls to which they are attached and, therefore, must be firmly secured in order to support the loads for which they are intended.

### Available Models

<b>RBA4007-312</b>	Drop Down Rail for Bariatric, Left Hand
<b>RBA4007-313</b>	Drop Down Rail for Bariatric, including Toilet Roll Holder, Left Hand
<b>RBA4007-314</b>	Drop Down Rail for Bariatric, Right Hand
<b>RBA4007-315</b>	Drop Down Rail for Bariatric, including Toilet Roll Holder, Right Hand
<b>RBA4007-316</b>	Drop Down Rail for Bariatric, including Nurse Call Button (Coms by Others), Left Hand
<b>RBA4007-317</b>	Drop Down Rail for Bariatric, including Nurse Call Button (Coms by Others), Right Hand

### Related Products

<b>RBA8851-137</b>	Bariatric Toilet
<b>RBA8186-853</b>	Big John Bariatric Toilet Seat, Closed Front

**Important:** Installation instructions and current rough-in details should be furnished with each fixture. Do not rough in without certified dimensions. Dimensions are subject to manufacturer's tolerance of +/-10mm.

